



A **covetrus**  Company



Modular Plating Box

Portability, Presentation & Personalisation



A **covetrus** Company

Veterinary Instrumentation Modular Plating Box

Portability, Presentation & Personalisation

Advancing the supply of Veterinary Instrumentation's most clinically current fracture fixation plating kits, we introduce the Modular Plating Box concept.

Starting with a universal basket and lid, this purpose built cassette can house an array of bespoke instrumentation trays, plate trays and screw caddies, all designed with tailor-made holding devices to accommodate plating kits in size ranges from 1.5mm to 4.5mm.

Portable in nature, the handling characteristics deliver a seamless processing experience. The concept ensures structure and formation, improving product identification and **presentation**.

With each plating box being modular in nature, the user has the ability to **personalise** and tailor the components adopted, offering users the opportunity to adopt elements in a supplementary fashion, or as part of a full plating kit.

Manufactured in the US, the concept is dynamic yet robust and provides protection to the components inside, ensuring longevity and quality are personified.

Compatible with all types of sterile packaging technologies, users can deploy disposable wrap or rigid sterile packaging and for components used in a supplementary manner, sterilisation peel pouches.

Underpinning the above is the efficacy of our implants, the crucial components of any plating kit.

With over **19,000** Vi products being implanted each year from a diverse portfolio of over **1,400** different plate and screw profiles, validation and user assurance is achieved through a revision rate of under **0.1%**, rivalling any fracture repair system.

 [@VeterinaryInstrumentation](#)

 [@Veterinary Instrumentation, a Covetrus Company](#)

T: 1800 785 642

E: vi@provet.com.au

W: www.veterinary-instrumentation.com.au

Veterinary Instrumentation, 48 Bell-Are Avenue,
Northgate, QLD 4013




A covetrus Company

Modular Plating Box Configurations:

Product Code & Description		
MPBL – Modular Plating Box Lid	305x130x10mm	0.22kg
MPBB – Modular Plating Box Base	305x130x115mm	0.56kg
<i>Combined</i>	315x145x120mm	0.78kg



Product Code & Description		
MPBIT1520 – 1.5/2mm Instrument Tray	230x130x35mm	0.29kg
MPBIT24 – 2.4mm Instrument Tray	230x130x35mm	0.31kg
MPBIT27351 – 2.7/3.5mm Instrument Tray 1	230x130x35mm	0.31kg
MPBIT27352 – 2.7/3.5mm Instrument Tray 2	230x130x35mm	0.27kg
MPBIT45 – 4.5mm Combined Instrument & Plate Tray	300x125x40mm	0.37kg



 [@VeterinaryInstrumentation](#)

 [@Veterinary Instrumentation, a Covetrus Company](#)

T: 1800 785 642

E: vi@provet.com.au

W: www.veterinary-instrumentation.com.au

Veterinary Instrumentation, 48 Bell-Are Avenue,
Northgate, QLD 4013



A **covetrus** Company

Product Code & Description		
MPBSC15 – 1.5mm Screw Caddy	70x60x30mm	0.076kg
MPBSC20 – 2mm Screw Caddy	70x60x42mm	0.082kg
MPBSC24 – 2.4mm Screw Caddy	70x60x42mm	0.082kg
MPBSC27 – 2.7mm Screw Caddy	70x60x42mm	0.082kg
MPBSC35 – 3.5mm Screw Caddy	128x76x65mm	0.184kg
MPBSC45 – 4.5mm Screw Caddy	128x76x65mm	0.184kg



For a detailed reference on what screw sizes and the volume of each which can be stored in each respective screw caddy, please refer to the supporting document by scanning the QR code below:



 [@VeterinaryInstrumentation](#)

 [@Veterinary Instrumentation, a Covetrus Company](#)

Veterinary Instrumentation, 48 Bell-Are Avenue,
Northgate, QLD 4013

T: 1800 785 642

E: vi@provet.com.au

W: www.veterinary-instrumentation.com.au



A covetrus Company

Product Code & Description		
MPBPTL15 – 1.5mm Plate Tray For Locking Plates	128x45x16mm	0.1kg
MPBPTNL15 – 1.5mm Plate Tray For Non Locking Plates	128x45x16mm	0.1kg
MPBPTL2 – 2mm Plate Tray For Locking Plates	128x45x16mm	0.1kg
MPBPTNL2 – 2mm Plate Tray For Non Locking Plates	128x45x16mm	0.1kg
MPBPTL24 – 2.4mm Plate Tray For Locking Plates	128x45x16mm	0.1kg
MPBPTNL24 – 2.4mm Plate Tray For Non Locking Plates	128x45x16mm	0.1kg
MPBPTL2735 – 2.7/3.5mm Combined Plate Tray	230x130x35mm	0.31kg



If you would like to find out more about the new Modular Plating Box range please contact Vi or your local market representative.



 [@VeterinaryInstrumentation](#)

 [@Veterinary Instrumentation, a Covetrus Company](#)

T: 1800 785 642

E: vi@provet.com.au

W: www.veterinary-instrumentation.com.au

Veterinary Instrumentation, 48 Bell-Are Avenue,
Northgate, QLD 4013



A **covetrus**  Company

