

## 1.5mm & 2.0mm Combi T-Plate

The Vi Combi T-Plates have been created as an extension of our current T-Plate range. The general shape of the plate is very similar, and this plate comes with the added **benefit of locking technology**. There are stacked locking holes in the distal end of the plate and dynamic locking holes in the shaft.

### Key Features & Benefits:

- Locking construct **increases stability**
- Stacked locking hole allows placement of either a locking or non-locking screw
- Dynamic Locking holes in proximal section
- K-wire hole
- Available in a range of lengths and 2 different sizes
- Available in Stainless Steel



### Product Ordering Details

Product Code	Product Description
T151532COMBI	1.5mm Combi T-Plate 32.5mm Long 6 Holes
T1515735COMBI	1.5mm Combi T-Plate 35.5mm Long 7 Holes
T202033COMBI	2.0mm Combi T-Plate 33mm Long 5 Holes
T202040COMBI	2.0mm Combi T-Plate 40mm Long 6 Holes
T2020743COMBI	2.0mm Combi T-Plate 43.5mm Long 7 Holes
T2020850COMBI	2.0mm Combi T-Plate 50.5mm Long 8 Holes